Neurociencia cognitiva

2do cuatrimestre 2024

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


Psychopy

1. Builder

The screenshot displays the PsychoPy Builder interface for a routine named "prueba.psyexp". The window title is "prueba.psyexp - PsychoPy Builder (v2024.2.1)". The interface includes a menu bar (File, Edit, View, Tools, Experiment, Demos, Pavlovia.org, Window, Help) and a toolbar with icons for File, Edit, Experiment, Desktop, Browser, and Pavlovia. The main workspace shows a timeline for the routine "prueba" with a time axis labeled "t (sec)" ranging from 0 to 3.5. Two components are visible: "prueba" (Text) and "key_resp" (Keyboard). The "prueba" component is active from t=1.0 to t=3.5, and the "key_resp" component is active from t=0.0 to t=3.5. A "Routine settings" panel is visible on the left. On the right, a "Components" panel shows a list of available components: Text, Keyboard, Sound, Mouse, Image, Slider, Stimuli, Responses, Custom, and EEG. The "Flow" section at the bottom shows a diagram with a "trial" box and a "trials" box, with an "Insert Loop" button highlighted.

Psychopy

2. Coder

The screenshot displays the PsychoPy Coder (v2024.2.1) interface. The main window shows a Python script named `atencion_probe.py` with the following code:

```
19 # Import libraries
20 from psychopy import visual, event, core, data, gui
21
22
23 # Genera la ventana en unidad pixels
24 win = visual.Window([1600, 900], fullscr=True, units='pix')
25
26 info = {} # a dictionary
27 # present dialog to collect info
28 info['participant'] = 'prueba_1'
29 dlg = gui.DlgFromDict(info) # (and from psychopy import gui at top of script)
30 if not dlg.OK:
31     core.quit()
```

The interface includes a menu bar (File, Edit, View, Tools, Demos, Pavlovia.org, Window, Help), a toolbar with icons for file operations and running the experiment, and a status bar at the bottom showing "Line: 1 Col: 1" and "Python". A shell window at the bottom left contains the prompt "Hit [Return] to start a Python session."

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